



Product Code . DS-LE-11619

Hand Lenses (Magnifying Glasses)

Description

Catalogue No. : “IDAC229/036”

DESCRIPTION :-

A magnifying glass (called a hand lens in laboratory contexts) is a convex lens that is used to produce a magnified image of an object.

The lens is usually mounted in a frame with a handle

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/