

Product Code . DS-LE-11657

Microscope (Light)



Description

Catalogue No. : "IDAC131/074"

Features :-

Precisely designed

Low maintenance

Hassle free functionality

Easily adjustable lenses

Replaceable components

Easy to carry

Excellent performance

Long service life

DESCRIPTION :-

Optical system: Infinite optical system

Viewing head: Seidentopf Binocular head, Inclined at 30 degree

Interpupillary Distance 55-75 mm

Eyepiece: WF10X/20

Infinite Semi-Plan Achromatic Objective

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/