



Product Code . DS-LE-11758

Gas Refrigerator

Description

Catalogue No. : "IDAC310/0175"

DESCRIPTION :-

Medical Refrigerator is specialized refrigerating equipment for cold

Foamed -in-place CFC free PUF insulation ensure temperature stability and reduced energy temperature.

Storage of pharmaceutical as well as logical products and so forth

Micro-compressor Control & Large screen LED Display

Temperature range: 2 - 8°C x 0.1°C

Audible and visual alarm for: High and low temperature, System & sensor failure, Door ajar, Power failure alarm, Low backup battery alarm.

Adjustable shelves made of high-quality steel wire, 5pcs shelves supplied

Inside door operated illumination is provided for visibility of bags kept inside the chamber.

Power Supply: AC220V±10%, 50/60Hz, 200W

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/