

Product Code . DS-LE-11767

Industrial Oven



Description

Catalogue No. : “IDAC111/0184”

DESCRIPTION :-

Double walled construction inner chamber made of Stainless Steel & exterior made of G.I. sheets powder coated.

The Gap between the two walls filled with glass wool insulation.

Fully insulated door with brass lock hinges & aluminum ventilator.

Heating is by specially designed heaters made of Nichrome wire fitted at the bottom of the Oven.

Temperature Controlled by Electronic Digital Temperature controller with PT 100 as sensor.

Supplied with G.I. wire mesh shelves.



Temperature Range: 50oC to 250oC

Temperature Resolution: 0.1oC

Control Accuracy: +-2oC

Material: Stainless Steel

Power Consumption: 220 V

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/