

Product Code . DS-CW-11773

Weight Box



Description

Features :-

- Brand: Didac Scientific
- Model: DAC146
- Country of Origin: United Kingdom
- OIML Class: E2
- Material: Stainless Steel
- Weight: 2110 g
- Nominal value: 1 g - 1 kg
- Quantity: 13 pcs

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/