

Description

Lab Bench Oven Manufacturers UK

The UK is a leader in manufacturing high-quality lab bench ovens Lab Bench Oven Manufacturers UK. These trusted companies offer a diverse range of ovens designed for various laboratory applications. From compact benchtop models for routine drying and sterilization to larger capacity ovens for demanding research projects, UK manufacturers cater to your specific needs. Explore their innovative features like precise temperature control, digital displays, and safety features for reliable and efficient performance in your laboratory.

Description:

Use: General drying of samples in the laboratory.

Interior body is made of stainless steel by argon-arc-welding technique.

Exterior body is made of high-quality steel sheet with a novel appearance.

Oven door with Silicon gasket & transparent thermal glass window for watching the samples under drying.

Temperature control protection.



Hot air circulation system is composed
of the blower fan and proper air
channels.

Digital display.

Microcomputer PID temperature
controller with a timer function to ensure
precise & reliable temperature control.

Can adjust blast and exhaust in control
panel.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108
Martinfield, Welwyn Garden City ,
United Kingdom

Direct Contact Details 
sales@didacscientific.co.uk
 www.didacscientific.co.uk/